



# Glass vessel


GF1000K with inlet pipe for SRK 270/160 probe



- Feed vessel for filling machines
- Sterilizable
- Level measurement with inlet pipe for impedance SRK bar probe
- Schott Duran 50 glass
- Total vol. 1000 ml / working vol. 637 – max. 774 ml

## Use

The sterilizable receiving vessel GF1000K for filling systems in the sterile area is operated with the level control system LLCU (mipromex MLT). The glass vessel is attached to the system with an Aquasant sheet steel holder or attached to the stand rod. This vessel ensures that the filling machine is constantly supplied with product.



## Technical Data

### Design

glass cylinder closed with tube inlet for SRK bar probe, tube connections for ventilation, inlet and outlet

### Installation

into aquasant® POM support arm art.no.: 04.40.33.010

### Function

glass receiver with filling level probe for continuous level monitoring, lateral inlet, central outlet connections

### Handling

connect respective silicone tubing and secure hose with clamp.

V1 safety overfilling volume 250 ml

V2 working volume 637 ml

V3 residual volume 137 ml

Weight 570 g

Material Schott Duran 50 glass

### Connection

probe SRK 270/160 SB R N GF art.no.: 02.29.12.0058  
MTI 100/. AG / Liquid level control unit LLCU, MLT6130

### Exhaust connection

ventilation connection 2 x ø 6/10

### Silicone tubing

ø10/15 ±0.5 mm, shore: 55° ±5

### Working pressure

unpressurized, max. 1 bar

### Use in sterile area A

### Cleaning temperature

210 °C, max 10 min unpressurized (CIP/SIP) sterilizable and autoclavable

### Label

product number, serial number and logo

### Packaging

wrapped in foil, not sterile

## Function

The impedance changes in function of the dielectric and electrically conductive properties of organic products and aqueous solutions, as well as the immersion depth of the active measuring electrode. The measured impedance is converted directly into a standardized digital signal by the MTI measuring electronics as a sum signal and transmitted as a pulse packet to the mipromex® MLT6130.

### Conformity

Conformity evaluation procedure according to module A category I

CE conformity to EN 62061 category: 1 / EN ISO 13849 PL: a

ISO 9001:2015 CE1254

